

# excellence for life



Winter  
Arrhythmia  
School

Blue Mountain

February 6-8, 2015

BIOTRONIK // WAS 2015: MRI in CIED

# **Myths and Facts: A Discussion of MRI Imaging for Pacemaker and ICD Patients**



**Matt Reitzel**

**BIOTRONIK**  
**Canada**

## Topics of discussion:

---

- **Your experience with MR-conditional devices**
- **Overview of MR-conditional devices in the field**
- **Components of MR-conditional systems**
- **CHRS recommendations for device/patient management for MR scan**

## Topics of discussion:

---

- **Your experience with MR-conditional devices**
- **Overview of MR-conditional devices in the field**
- **Components of MR-conditional systems**
- **CHRS recommendations for device/patient management for MR scan**

## **Clicker Question: Does your centre currently implant MR-conditional CIED?**

---

- A) No, 0%**
- B) Yes, <5%**
- C) Yes, <20%**
- D) Yes, > 20%**

**Clicker Question: Does your centre currently provide on-site MR scans for patients implanted with MR-conditional CIED?**

---

- A) No scanning, no protocol to scan**
- B) No scanning, protocol in place to scan if required**
- C) Yes scanning preformed and protocol in place**

## Topics of discussion:

---

- Your experience with MR-conditional devices
- **Overview of MR-conditional devices in the field**
- Components of MR-conditional systems
- CHRS recommendations for device/patient management for MR scan

# BIOTRONIK ProMRI® - Developing the industry's most complete MR-conditional CRM portfolio

**BIOTRONIK**  
starts  
development of  
MR-Conditional  
systems



**Approval of 1<sup>st</sup>  
MR Cond. ICD  
System**



**Launch MR Cond.  
ICD System**



**BioMonitor**



**Expansion MR  
Cond. ICD  
System**



**MR Cond. FBS &  
3.0 T ICD System**



**BIOTRONIK's  
1<sup>st</sup> MR Cond.  
IPG System**



**Expansion of  
the MR Cond.  
IPG Portfolio**



**MR Cond. CRT-  
P Systems**



**MR Cond. Full  
Body Scan IPG  
System**



**Next generation  
MR Cond. IPG  
System**

**2003**

**2008**

**2009**

**2010**

**2011**

**2012**

**2013**

**2014**

**Approval of 1<sup>st</sup> MR  
Cond. IPG System**

Document identifier

**Launch of 1<sup>st</sup> MR  
Cond. IPG System**



# Topics of discussion:

---

- Your experience with MR-conditional devices
- Overview of MR-conditional devices in the field
- **Components of MR-conditional systems**
  - What's the difference? Combination of leads/devices? Exclusion zone or full body scan?
- CHRS recommendations for device/patient management for MR scan
- Future of MR-conditional devices

**Clicker Question: It is estimated that a patient has an X % chance of requiring an MR examination over his/her lifetime after CIED implantation. X = ?**

---

- A) 0-25%**
- B) 26-50%**
- C) 51-75%**
- D) 75-100%**

**Clicker Question: It is estimated that a patient has an X % chance of requiring an MR examination over his/her lifetime after CIED implantation. X = ?**

---

**A) 0-25%**

**B) 26-50%**


**C) 51-75%**

**D) 75-100%**

Kalin R, Stanton MS. Current clinical issues for MR scanning of pacemaker and defibrillator patients. Pacing Clin Electrophysiol 2005;28:326-8.

## 50-75% of patients with a CIED are indicated for an MR scan


- Over 30 million MRI scans conducted in the US<sup>1</sup>.
- 2.5 million US patients have a pacemaker, and this number is growing by 80.000 per year<sup>2</sup>.
- Estimated 50%-75% of patients with a device are indicated for an MRI scan<sup>3</sup>.



# Utilization of MRI Scan

## MRI utilization increases by age

- Patients aged 70-74 years are most likely to receive an MRI scan\*



\* BARMER Health insurance: Arztreport 2011

# Do the terms MRI Safe and MRI conditional have the same meaning?

---



**MR- safe:** An item that poses no known hazards in all MRI environments



**MR- conditional:** An item that has been demonstrated to pose no known hazards in a specified MR environment with specified conditions of use (field strength, spatial gradient etc.)



**MR-unsafe:** An item that is known to pose hazards in all MR environments

# MRI compatibility- what's the big deal?

---



# Manufacturing MRI compatible devices- What's the difference?

- Three Birdcages for 1.5 & 3.0 Tesla testing
- Shielded chambers
- RF antenna of clinical MRT is used



16

Document identifier



# Manufacturing MRI compatible devices- What's the difference? Software

- Spec
- Asyn
- durin
- system
- Non-
- depe
- Diag
- meas
- arrhy
- oper

Parameters | Iilesto 7 DR-T XXXXXXXX

ICD therapy  
OFF ON

Prog.

Follow-up

Parameters

Tests

Recordings

Diagnostics

Status

More

Preferences

End

### MRI checklist

**Check device and leads**

- A dedicated BIOTRONIK MR conditional lead and device form an MRI-tested system
- Device has been implanted in the pectoral region for more than 6 weeks
- Follow-up was successful and threshold does not exceed 2.0 V at 0.4 ms
- No additional implanted devices present
- No additional leads, wire adapters or lead extenders present

**Radiological considerations**

- Standard 1.5 T cylindrical scanner architecture required
- Continuous patient monitoring required during MR scan
- Observe specific conditions for MR Conditional devices (SAR, scan zone, field strength, slew rate ...)
- After MRI scan, restore previously programmed parameters, program and confirm settings.

MRI program ON

MRI mode OFF

Basic rate [bpm]

☒ **accept the conditions for MRI examinations**

- Please select an MRI mode. After clicking 'OK' the MRI parameters will be displayed.
- Selecting 'Program' on the next screen will program the device accordingly.
- Any parameter change will result in the loss of the MRI program.

Print Help OK Cancel

Supporting

ing

# **Manufacturing MRI compatible devices- What's the difference? Hardware**

---

- **Improved Electromagnetic Interference (EMI) susceptibility**
- **New chip developed to improve EMI rejection and minimize the number of external component**
- **Optimized component selection and layout for improved EMI rejection over 16Hz - 3000Mhz**
- **Significantly reduced the length of the electrical connections between the feedthrough and the electronic circuit**
- **Minimized EMI radiation inside pacemaker housing**

**Clicker Question: Approved ProMRI systems are MR-safe.**

---

**A) True**

**B) False**

## Clicker Question: Approved ProMRI systems are MR-safe.

---

**A) True**

**B) False**

No CIED system is "MR-safe", but selected CIED systems are "MR-conditional"

**Clicker Question: Conditions for MR imaging of an MR-conditional CIED are consistent from vendor to vendor and device to device.**

---

**A) True**

**B) False**

## **Clicker Question: Conditions for MR imaging of an MR-conditional CIED are consistent from vendor to vendor and device to device.**

---

**A) True**

**B) False**

Manufacturer-specified scanning parameters are always checked and followed; the MR scanning protocol will vary in accordance with the patient's CIED system. Some systems are initially approved with exclusion zone and then retroactively approved for full body scanning. Recommended to check with vendor before scan to ensure parameters are known.

## **Clicker Question: MRI for a device with an exclusion zone.**

---

**The manufacturer guidelines state that the MRI isocentre can pass no lower than eye level and no higher than hip level. Are images of the C and L spine possible?**

**A) Yes**

**B) No**

## Clicker Question: MRI for a device with an exclusion zone

---

**Are images of the C and L spine possible?**


**A) Yes**

**B) No**



# MR Exclusion Zone


## Exclusion zone in daily practice



# MR Exclusion Zone

## Exclusion zone in daily practice

Scan with isocenter on eye level




isocenter eye level

# MR Exclusion Zone

## Exclusion zone in daily practice

Scan with isocenter on eye level



isocenter eye level

**C2  
C3  
C4  
C5  
C6  
C7  
T1  
T2  
T3  
T4  
T5**


**Head, cervical spine and upper thoracic spine (T4) can be examined**

**Real patient case with isocenter on eye level and Siemens MRI scanner. Patient length 1.75 m and scan taken within the ProMRI conditions.**

# MR Exclusion Zone

## Exclusion zone in daily practice

Scan with iso-center on hip level



**Lower lumbar spine can be examined**

**Real patient case with isocenter on hip level and Siemens MRI scanner. Patient length 1.75 m and scan taken within the ProMRI conditions.**

\* depends by MRI model

## **Clicker Question: A device with an MRI exclusion zone limits what percentage of total indicated MRI Scans?**

---

- A) <5%**
- B) 10-15%**
- C) 15-25%**
- D) 25-50%**
- E) >50%**

## Clicker Question: A device with an MRI exclusion zone (thoracic area) limits what percentage of total indicated MRI Scans?

---

A) <5%

**B) 10-15%**

C) 15-25%

D) 25-50%

E) >50%


BARMAR: Arztreport 2011. Barmar GEK- large Health Care provider data in Europe

# 88% of MRI scans can be performed

- BARMER health insurance data show the MRI usage per scan location
- Data are verified by BIOTRONIK radiology consultants

Head	26,3%
Neck	2.0%
Cervical & Lumbar spine	26.7%
Hip	3.2%
Extremities	30.2%

**88.4 %**




\*BARMAR: Arztreport 2011

\*\* C, T and L spine percentages based on expert evaluation

**Clicker Question: Are radiopaque markers on devices and leads a reliable way to determine MR-compatibility?**

---



**A) Yes**

**B) No**



## **Clicker Question: Are radiopaque markers on devices and leads a reliable way to determine MR-compatibility?**

---

**A) Yes**

**B) No**

How would we tell only with a marker if it was full body, exclusion zone, 1.5T or 3T etc. Also, as “legacy” leads are approved for MR-compatibility, these markers would not appear on these leads. It is unreliable due to this case and would also require checking with the manufacturer for the acceptable scanning parameters when these markers are present

## **Clicker Question: Biotronik has “back dated” over 30 non-ProMRI lead types as MR Conditional (CE Mark)**

---

**With CE Mark, if a ProMRI pacemaker is implanted with these leads, is the system MRI compatible in Canada?**

**A) Yes**

**B) No**

## **Clicker Question: Biotronik has “back dated” over 30 non-ProMRI lead types as MR Conditional (CE Mark)**

---

**With CE Mark, if a ProMRI pacemaker is implanted with these leads, is the system MRI compatible in Canada?**

**A) Yes**

**B) No**

Along with the current leads, BIOTRONIK has tested 34 current, non-ProMRI leads and approved them as MRI conditional (with CE Mark). At this time, Health Canada has not granted the back date and therefore, only the currently approved system is considered ProMRI

**Clicker Question: A patient is implanted with MR-compatible leads from Company X and an MR-compatible pacemaker from Company Y.**

---

**Is this system MR-compatible?**

**A) Yes**

**B) No**

**Clicker Question: A patient is implanted with MR-compatible leads from Company X and an MR-compatible pacemaker from Company Y.**

---

**Is this system MR-compatible?**

**A) Yes**

**B) No**

Labeling for Health Canada has approved "device systems". Although individual components are MR-compatible, this has not been validated and stated as an MR-conditional system.

## Topics of discussion:

---

- Your experience with MR- conditional devices
- Overview of MR- conditional devices in the field
- Components of MR- conditional systems
- **CHRS recommendations for device/patient management for MR scan**

## **Clicker Question: Which of the following are absolute contraindications for MR scanning?**

---

- A) MR scanning protocol outside of manufacturer specifications**
- B) Presence of abandoned pacing or ICD leads**
- C) Broken or fractured lead(s)- known or suspected**
- D) Permanent epicardial pacing or ICD leads**
- E) A,B and C**
- F) All of the above**

## **Clicker Question: Which of the following are absolute contraindications for MR scanning?**

---

- A) MR scanning protocol outside of manufacturer specifications**
- B) Presence of abandoned pacing or ICD leads**
- C) Broken or fractured lead(s)- known or suspected**
- D) Permanent epicardial pacing or ICD leads**
- E) A,B and C**
- F) All of the above**



**Clicker Question: Which of the following methods is acceptable to monitor the patients vitals while undergoing MR scan?**

---

- A) pulse oximetry**
- B) Electrocardiographic (ECG)**
- C) Capability for verbal communication between the MR scan operator and patient**
- D) All of the above**

**Clicker Question: Which of the following methods is acceptable to monitor the patients vitals while undergoing MR scan?**

---

- A) pulse oximetry**
- B) Electrocardiographic (ECG)**
- C) Capability for verbal communication between the MR scan operator and patient**
- D) All of the above**

It is recommended to use either pulse oximetry or ECG in conjunction with verbal interaction according to CHRS consensus

## **Clicker Question: Monitoring of the MR-conditional CIED patient during MR scanning may be preformed by:**

---

- A) MR Technologist**
- B) MR Radiologist**
- C) CIED Clinic nurse**
- D) Cardiologist with knowledge of CIED**
- E) All of the above**

## **Clicker Question: Monitoring of the MR-conditional CIED patient during MR scanning may be preformed by:**

---

**A) MR Technologist**

**B) MR Radiologist**

**C) CIED Clinic nurse**

**D) Cardiologist with knowledge of CIED**

**E) All of the above**

According to CHRS consensus, all of the above are acceptable however, should be evaluated based on the patient's medical status, functional status of the CIED and experience of the staff

**Clicker Question: Current approved devices require a CIED clinic staff to test and reprogram devices pre- and post-MR scan.**

---

**A) True**

**B) False**

**Clicker Question: Current approved devices require a CIED clinic staff to test and reprogram devices pre- and post-MR scan.**

---

**A) True**

**B) False**

CIED clinic staff ensures device testing is within appropriate ranges for MR-scan (lead impedance, sensing and capture) and ensures MRI programming for the device during scan. Post scan testing to ensure device functionality and restores permanent programming.

## **Clicker Question: Next generation MR devices do not required a CIED clinic staff to test and reprogram devices pre- and post-MR scan?**

---

Note: Next generation of MR compatible devices may include devices that auto- detect MRI environment and reprogram themselves or have an MRI activator the radiologist may use to temporarily reprogram the device pre- and post- MR scan.

**A) True**

**B) False**

**Clicker Question: Next generation MR devices do not required a CIED clinic staff to test and reprogram devices pre- and post-MR scan?**

---

**A) True**

**B) False**

Regardless of MRI programming of the device, an full interrogation and evaluation of the CIED system is required pre- and post-MRI as per CHRS guidelines.



# Thank you for your participation!

---

